



**PRODUCT CATALOG**



## Innovative solutions in urology, skin protection and securement

BioDerm® is one of the leading manufacturers of proprietary hydrocolloid products in the United States for skin protection, urinary management, and securement.

The company challenges accepted inferior standards of care by creating products that reduce infection rates, add comfort and reliability, and vastly improve quality of life.

### UROLOGY

Men's Liberty™ and Men's Liberty Acute™ are revolutionary alternatives to diapers, condom catheters and when medically unnecessary, Foley catheters.




### SKIN PROTECTION

Ensure proper product application and removal with the use of BioPlus+™ Barrier Film Wipes and FreeDerm® Adhesive Remover.



### SECUREMENT

The CathGrip® family of universal medical tubing securement options provides clinicians with the ability to manage a wide variety of catheters and tubes ranging in size from 6-42 Fr.







# UROLOGY

## Say Goodbye To Adult Diapers Once And For All

Staying dry and odor-free shouldn't even be a question. The Men's Liberty line of products are non-invasive, external urine collection devices that can be used at home and on-the-go, so you can forget about bladder leakage and get back to living life on your own terms.

- Replaces diapers, condom catheters and when medically unnecessary, Foley catheters
- Has a completely external design that fits most male anatomy
- Keeps men dry and comfortable for up to 24 hours
- 8 oz collection pouch can be emptied easily throughout the day
- Connects directly to a leg or bed bag
- Covered by most insurance policies in the US (Standard deductibles apply)




### How Does It Work?

Men's Liberty & Men's Liberty Acute are made from hydrocolloid, an adhesive that creates a secure, skin friendly seal around the urinary opening. Men's Liberty is designed for use in the outpatient setting and includes an integral collection chamber that holds up to 8 fl. oz (250 mL) of urine. Patients can fill the collection chamber, which can then be discreetly emptied into a toilet through the end plug located on the bottom of the product, whenever they get the opportunity.

#### Position & Apply




#### Secure & Seal





## For Outpatient Use



### Men's Liberty™

Men's Liberty is the completely external catheter that features a micro-vent and is applied just once a day and keeps men dry and comfortable for up to 24 hours. This discreet, one-size-fits-most option manages men's urinary incontinence in a safe, healthy way. More than 5 million units have been sold. Men's Liberty is covered by Medicare and over 3,000 insurance plans in the US.

## For Acute Care Use



### Men's Liberty™ Acute


Men's Liberty Acute is a male external catheter designed for institutional and hospital use. It is completely external and seals gently and securely to the tip of the penis with BioDerm's proprietary hydrocolloid adhesive.



### Penile Clamp for Men

The KindKlamp, penile clamp, is designed to temporarily stop the flow of urine through the urethra when applying Men's Liberty. This simple accessory helps men who are dribbling or leaking and are having trouble applying the Men's Liberty device.

- Easy use clasp – can be closed and released with one hand!
- Curved upper arm for improved comfort
- Long lasting hinge designed for multiple cycles
- Designed to fit most male anatomy
- The plastic design allows for use during x-rays and disgnostic imaging



SKU# 20003

### KindKlamp™

#### For Outpatient Use

The KindKlamp is made from durable plastic and ultra soft foam for added comfort. It is made without wires or metal. This ensures that the KindKlamp is skin friendly and can be wiped with a damp cloth without the potential of rust or corrosion. Product use varies, but clamp should be replaced every 90 days.



**1000 mL Leg Bag**  
**SKU# 90000**

## Leg Bag

The Leg Bag can be used for extra capacity while using the Men's Liberty. The universal tube directly connects to most catheter systems. To use, simply unplug the bottom of the Men's Liberty pouch and connect to the Leg Bag.

- Kink-resistant tubing
- Convenient drain spout
- Anti-reflux valve
- All fluid pathways are sterile
- Straps included

*(US only)*



**2000 mL Bed Bag**  
**SKU# 90001**

## Bed Bag


The Bed Bag can be used for extra capacity while using the Men's Liberty during the night. The universal tube directly connects to most catheter systems. To use, simply unplug the bottom of the Men's Liberty pouch and connect to the Bed Bag. Two bags are covered by insurance every month with a doctor's order.

- Kink-resistant tubing
- Convenient drain spout
- Anti-reflux valve
- All fluid pathways are sterile
- Hook included to connect to bed railing

*(US only)*



# SKIN PROTECTION



## Faceplate Strip

Faceplate Strip prevents skin injuries and irritation due to face masks, glasses, and face shields. Face masks, glasses and face shields are a necessity and are being worn more frequently and for longer periods of time. No matter how well your face PPE fits, constant wear time can cause skin injuries and irritation. The wing design will cover the bridge of your nose, cheek bones and forehead for up to 24 hours.

- Made in the USA and are skin friendly
- Can be worn for up to 24 hours



Box of 25, SKU# 52300

## Bio Plus+<sup>TM</sup> Barrier Film Wipe

BioPlus+ improves hydrocolloid adhesion by up to 50%. A BioPlus+ wipe is included in most BioDerm products and is also available as a stand alone product.

- No-sting
- Forms transparent, protective barrier film
- Dries faster than major competitors in lab tests
- Increases skin adhesion more than major competitors in lab tests



## Adhesive Remover

Specially formulated to instantly remove hydrocolloid, tape and bandages.

Instantly removes hydrocolloid without irritation and evaporates rapidly with no residue. Allows immediate re-application of an adhesive device without washing. Reduces skin trauma from tape and dressing removal including infection risk, inflammatory damage and edema. Available as 1 oz or 3 oz pump spray and single-use wipes.

- Safe and gentle for all patients
- No-sting
- Instantly removes hydrocolloid
- No harsh chemicals
- Oil and lotion free

**Dries in 15 seconds**  
allowing for immediate reapplication



1 oz. SKU# 52200

3 oz. SKU# 52202

**NEW!**



6 mL Wipe

SKU# 90037

# SECUREMENT

---



## The CathGrip Family of Universal Securement Options

The CathGrip line provides clinicians with the ability to manage a wide variety of catheters and tubes ranging in size from 6-42 Fr. Each product is made of BioDerm's proprietary hydrocolloid with features that solve a wide variety of challenges. When properly applied, these devices provide a comfortable fit and secure hold for up to seven days. The soft straps are durable for repeat positioning, perfect for continual ease of use. CathGrip can meet all your securement needs, making it the new standard of care sought after by physicians, nurses and purchasing teams.




**Not made with  
natural rubber latex**



**Up to 7-day  
wear time**



**Fits 6-42 Fr.  
medical tubing**








**The CathGrip family of products is exclusively  
distributed in the US by Medline Industries.  
Contact Medline at [customerservice@medline.com](mailto:customerservice@medline.com)  
or 1-800-MEDLINE.**



# CathGrip Single Strap & Double Strap

CathGrip is a truly universal system that securely holds tubes from 6-42 Fr. Available in single or double straps, CathGrip provides all over securement with a skin-friendly hydrocolloid adhesive to prevent skin injuries. Soft straps are durable for repeat positioning.



- Available in single or double straps in small, medium, and large sizes
- Hydrocolloid adhesive oval
- Not made with natural rubber latex
- Soft straps with no-slip grip
- Simple repeat positioning
- Fits medical tubing from 6-42 Fr
- Durable repeat positioning
- Extended wear time, 2-3 times longer than competitors

## Available in Sterile & Non-Sterile




# CathGrip Low Profile Single Strap & Double Strap

A new addition to the CathGrip family. The Low Profile CathGrip has all of the qualities of the existing Single and Double Strap but with the addition of low profile straps. The low profile straps hug the body reducing the risk of snagging for a more comfortable wear. Sterile and non-sterile available.



# CathGrip Oval

Designed to prevent pressure injuries and skin tearing with low-profile proprietary hydrocolloid securement, especially important when dealing with tender skin. The ovals are sold in pairs. The low profile design allows for application in difficult areas. Sterile and non-sterile available, except for PN 54002 only available as non-sterile.



Uses listed above are not all inclusive. May vary based on tube location, infection prevention protocols and presence of exudating wounds.







# CathGrip PEG Protect™

Designed specifically for use with PEG tubes, PEG Protect brings our skin friendly hydrocolloid securement to long and short term PEG patients. Only available as sterile.

- Allows for regular tube rotation for traditional and low profile PEG tubes
- Reduces the risk of dislodgement and tube migration
- Provides a skin friendly hydrocolloid barrier against bodily fluid, bruising and excoriation
- Reduces potential for tube kinking
- Two piece design allows for variable securement location related to the PEG tube insertion site

**Secures both standard and low-profile PEG tubes: allowing regular tube rotation for traditional and low-profile PEG tubes, Nephrostomy Chest Tube, PEG Tube and J Tube**




SKU# 51501

# CathGrip Extra Small


Extends the security of our hydrocolloid base to support insertion site securement. CathGrip Extra Small securely holds any tube from 3-24 Fr. and can help reduce dislodgement and minimize the risk of infection.

- Universal securement for catheters, drains & tubes 3-24 Fr.
- Durable soft strap with no-slip grip for repeat positioning
- Skin friendly hydrocolloid material moves with the body
- Can be applied to the face & body
- Reduces the risk of tube dislodgement
- Easy application and removal

**Available in Sterile & Non-Sterile**



Uses listed above are not all inclusive. May vary based on tube location, infection prevention protocols and presence of exuding wounds.



# Now That You Know Our Products

BioDerm provides innovative solutions for customers around the world. Safer, easy to use and cost effective, our proprietary hydrocolloid, skin prep, adhesive remover and soft strap technology is the next evolution in patient care. Join the future of healthcare technology today!

## OUTPATIENT

BioDerm's products are sold through our customer service center, based in Largo, Florida. Our friendly customer care team is trained to understand the conditions associated with urinary incontinence and are waiting for your call.

## ACUTE CARE

BioDerm's CathGrip and Men's Liberty Acute products are sold exclusively through Medline who targets acute care, long-term care and other surgical centers in the U.S. The CathGrip line of products is transforming care in critical care units, interventional radiology, medical/surgical units, wound care clinics and more with solutions that help reduce the incidence of medical device related pressure injuries and catheters becoming dislodged.

## HEALTHCARE PROVIDERS

BioDerm works seamlessly with many HCP offices to help patients obtain insurance coverage for Men's Liberty. Contact us today for a HCP Provider Kit.

## INTERNATIONAL

BioDerm's products are currently available for sale in over 10 countries through a network of independent distribution partners in each country.

## PARTNERSHIPS

BioDerm products are available for private label collaboration for target markets and partnerships. Contact our care team today.



# ADDITIONAL INFORMATION

## Regulatory Information

BioDerm is registered with the US Food and Drug Administration and other international regulatory bodies. Copies of registration documents are available upon request or online at <https://bioderminc.com/library/#quality-documents>

Product Category	Role
Tape and bandage, adhesive – BioPlus+ Barrier Film Wipe	Specification Developer
Device, intravascular catheter securement – CathGrip Double Strap; CathGrip Oval; CathGrip Single Strap; CathGrip PEG Protect; CathGrip Extra Small; Tube Securement Device; Low Profile CathGrip Single Strap; Low Profile CathGrip Double Strap; Medium LVAD Securement Device; Large LVAD Securement Device	Contract Manufacturer; Manufacturer
Solvent, adhesive tape – FreeDerm Adhesive Remover; FreeDerm Bottles; FreeDerm Wipes	Specification Developer
Device, paste-on for incontinence, non-sterile - Faceplate Seals; Men's Liberty Acute with CathGrip; Men's Liberty Acute with Cathgrip and Towel; Men's Liberty Acute without CathGrip; Men's Liberty; Men's Liberty Acute, BioDerm XLS (Oval)	Contract Manufacturer; Manufacturer
Clamp, penile - BioDerm Penis Clamp; KindKlamp	Specification Developer
Bandage, elastic - Faceplate Strip	Manufacturer



\*does not apply to leg bag, bed bag and sterile products.

### US Food and Drug Administration

Establishment Registration Number:  
1063299  
BioDerm, Inc., 12320 73rd Court North,  
Largo, FL 33773  
Status: Active

### Authorized Representative: Europe

Advena Ltd  
Tower Business Centre,  
2nd Flr, Tower Street,  
Swatar, BKR 4013 Malta  
[www.advenamedical.com](http://www.advenamedical.com)

### UK Responsible Person : United Kingdom

Advena Ltd  
Pure Office, Plato Close,  
Tachbrook Park Warwick.  
CV34 6WE. United Kingdom  
[www.advenamedical.com](http://www.advenamedical.com)

### British Standards Institute

Certificate of Registration:  
Quality Management  
System ISO-13485:2016  
Certificate #: FM 565097  
Originally Registered: 07/06/2011

Product Name	HCP Code	Catalog/SKU Number*	Size	Units Per Box	Boxes Per Case	Expiration Date (from the date of manufacture)*
<b>UROLOGY</b>						
Men's Liberty	A4326	23046	One Size	35	12	36months*
Men's Liberty Acute	NA	21026	One Size	35	12	36months*
KindKlamp	A4356	20003	One Size	12	n/a	n/a
Drainage Bag – Leg	A4358	90000	1000 mL	50	4	60 months
Drainage Bag – Bedside	A4357	90001	2000 mL	50	4	60 months
Faceplate Seals (Pack of 10)	NA	20002	One Size	1 Pack	n/a	48 months
<b>SKIN PROTECTION</b>						
BioPlus+ Wipe	A5120	52300	One Size	25	12	36 months
FreeDerm Spray (1oz)	A4455	52200	1 oz.	36	4	36 months
FreeDerm Spray (3oz)	A4455	52202	3 oz.	63	1	36 months
FreeDerm Wipe (6mL)	A4455	90037	6 mL	30	n/a	36 months
Faceplate Strip	NA	20002NP	4.5" x 1"	10	100	48 months
<b>SECUREMENT - CATHGRIP</b>						
Double Strap Large (2 Straps)	A4333 for Urinary	51300	4" x 4.75"	30	3	36 months*
Double Strap Large (2 Straps, Non-Sterile)		51300NS	4" x 4.75"	30	3	36 months*
Single Strap Large (1 Strap)		51350	4" x 4.75"	30	3	36 months*
Single Strap Large (1 Strap, Non-Sterile)		51350NS	4" x 4.75"	30	3	36 months*
Double Strap Medium (2 Straps)		51301	3" x 3.5"	30	3	36 months*
Double Strap Medium (2 Straps, Non-Sterile)		51301NS	3" x 3.5"	30	3	36 months*
Single Strap Small (1 Strap)		51250	3" x 3.5"	30	3	36 months*
Single Strap Small (1 Strap, Non-Sterile)		51250NS	3" x 3.5"	30	3	36 months*
Low Profile Single Strap Small (1 Strap)		51250LP	3" x 3.5"	30	3	36 months*
Low Profile Double Strap Large (2 Straps)		51300LP	4" x 4.75"	30	3	36 months*
Low Profile Double Strap Medium (2 Straps)		51301LP	3" x 3.5"	30	3	36 months*
PEG Protect	A5200 for Percutaneous	51501	3" x 3.5"	30	3	36 months*
Oval Large	Varies Based on Type of Tubing Secured	54003	3" x 3.5"	30	3	36 months*
Oval Large (Non-Sterile)		54003NS	3" x 3.5"	30	3	36 months*
Oval Medium		53002	1.48" x 2.4"	30	3	36 months*
Oval Medium (Non-Sterile)		53002NS	1.48" x 2.4"	30	3	36months*
Oval Small (Non-Sterile)		54002NS	.93" x 1.5"	30	3	36 months*
Extra Small		51023NW	1.48" X 2.4"	30	3	24 months*
Extra Small (Non-Sterile)		51023NSNW	1.48" X 2.4"	30	3	24 months*

\*Expiration date is based on the earliest expiration date of the components. CathGrip, Men's Liberty and Men's Liberty Acute ship without alcohol wipes in Canada, EU and UK. Product availability is dependent on the regulations of your country. Not all products are available in all countries.

A photograph of three medical professionals in a clinical setting. In the foreground, a woman with brown hair tied back, wearing a white lab coat and a blue stethoscope, looks directly at the camera with a slight smile. Behind her, two men in white lab coats are looking at a tablet. The man on the left has a beard and is gesturing with his hand. The man on the right is older, wearing glasses and a red tie. The background is a blurred hospital or clinic environment.

## Contact Us

### **For Our U.S. Customers**

Call 1-800-814-3174

Email [CustomerCare@BioDermInc.com](mailto:CustomerCare@BioDermInc.com)

Fax completed LMNs to (888) 920-9370

### **For International Inquiries**

Please contact [CustomerCare@BioDermInc.com](mailto:CustomerCare@BioDermInc.com)



12320 73rd Court North, Largo, FL 33773

US Sales: 800-814-3174

International Sales: +01 727-507-7655

Support: 800-864-6755

Email: [CustomerCare@BioDermInc.com](mailto:CustomerCare@BioDermInc.com)

Web: [BioDermInc.com](http://BioDermInc.com)

© 2023 BioDerm, Inc. All rights reserved. BioDerm®, FreeDerm® and CathGrip® are registered trademarks of BioDerm, Inc. BioPlus+™, KindKlamp™, PEG Protect™, Men's Liberty™ and Men's Liberty™ Acute are trademarks of BioDerm, Inc. Other trademarks acknowledged.  
Patent Info: <https://bioderminc.com/patents/>

LT 93943 H