

# ERADICATE MRSA IN AS LITTLE AS 4 HOURS\*

SILVERLON® LIFESAVER™  
IV CATHETER DISCS

SILVERLON® Lifesaver™ IV Catheter Discs are an antimicrobial dressing with a unique permanent silverplated metallic surface ideal for IVs, central lines, pin sites, VAD drivelines, ECMO Cannulas and chest tubes. Utilizing a large reservoir of silver ions, Lifesaver Discs continuously deliver an effective level of silver ions in the dressing for up to 7 days.



Silverlon® Lifesaver™  
IV Catheter Dressing  
1" 4.0mm I.D.

**SILVERLON®**


PERFORMANCE. POWER. SPEED.

Silverlon.com | 888.508.2872 or 630.232.2510

# SILVERLON® HAS AN ANTIMICROBIAL EFFECT WITHIN 2 HOURS AGAINST MRSA IN THE DRESSING, WITH ERADICATION IN 4 HOURS\*.

## IN VITRO TEST DATA

### SILVERLON® MRSA and PSEUDOMONAS KILL RATE



\*Citation: 510(k) K122817, page 174 of 379

BACTERIAL SPECIES	Time and CFU/ml							
	0.0 <sup>hr</sup>	0.5 <sup>hr</sup>	1.0 <sup>hr</sup>	2.0 <sup>hr</sup>	4.0 <sup>hr</sup>	6.0 <sup>hr</sup>	24.0 <sup>hr</sup>	
Staphylococcus aureus MRSA ATCC 33591	8,250,000	2,600,000	270,000	500	0	0	0	
Pseudomonas aeruginosa ATCC 9027	9,000,000	4,000,000	0	0	0	0	0	

## PERFORMANCE

Silverlon has been studied in multiple specialties, including an independent, peer-reviewed published study covering over 60,000 central line days. To view complete studies, visit [Silverlon.com](http://Silverlon.com).

## POWER

Lifesaver™ Catheter Discs are effective against a broad spectrum of pathogens and contain up to 100 times more silver than impregnated silver dressings.

## SPEED

SILVERLON® produces a broad spectrum antimicrobial effect that rapidly kills pathogens in the dressing, including MRSA and Pseudomonas, in as little as four hours.

### IVCD-MINI



Catheter  
5/8" Dressing  
1.5mm I.D.

### IVCD-15



Catheter  
1" Dressing  
1.5mm I.D.

### IVCD-40



Catheter  
1" Dressing  
4.0mm I.D.

### IVCD-70



Catheter  
1" Dressing  
7.0mm I.D.

### IVCD-15-85



Catheter  
1.5" Dressing  
8.5mm I.D.

### DPD-20-12



Catheter  
2" Dressing  
12mm I.D.

**SILVERLON®**

PERFORMANCE. POWER. SPEED.

[Silverlon.com](http://Silverlon.com) | 888.508.2872 or 630.232.2510